


# PEDO FORCEPS

This is only a sampling of the pedodontic extracting forceps that we carry. If you do not see the design you like, please ask, we will probably be able to get it for you. All of our extracting forceps are manufactured in Germany and have a lifetime warranty against defects in material and workmanship. They may be serrated at no charge at the time of purchase.




Universal, lower incisor & bicuspid  
**200-073** No. 151sk with spring  
**200-071** No. 151s no spring



Upper roots  
**200-084** No.51s with spring  
**200-095** No. 51s no spring



Lower cowhorn  
**200-075** No. 23sk with spring  
**200-076** No. 23s no spring



Universal, Incisor & bicuspid  
**200-086** No. 150as Upper, Bifurcated  
**200-087** No. 151as Lower, Bifurcated  
**200-085** No. 101as Bifurcated



Lower molars  
**200-074**  
 No. 17sk with spring.



Upper bicuspids  
**200-083**  
 No. 139 with spring



est. 1954

# E.A. Beck & Co. Tip Edge Course

Committed To Exceeding Your Expectations!

## PLIERS

- Highest quality surgical steel
- Lifetime warranty against defects
- Many patterns to choose from
- Finest craftsmanship

SEE ORDER SHEET FOR  
 SPECIAL PRICING FOR  
 COURSE PARTICIPANTS



### TWEED ARCH

To be used for positive control over wires to be held or torqued. Blade thickness is .045". Jaws parallel on edgewise wire. OTHER PATTERNS AVAILABLE.  
**Order Code: 500-006**



### STANDARD BIRD BEAK

An excellent wire bending plier for use on round or rectangular wire up to .030" in diameter. OTHER PATTERNS AVAILABLE  
**Order Code: 500-001**  
 Made in Poland  
**Order Code: 560-001**  
 Made in USA



### HOLLOW CHOP PLIER

Ideal for reshaping and contouring NiTi™ wire and arches and cuspid bands.  
**Order Code: 570-128**  
 Made in Japan  
**Order Code: 570-128p**  
 Made in Poland



### CINCHING PLIER

Small, serrated, tapered tips are ideal for hard to reach areas. Excellent for crimping or forming arch wires up to .028".  
**Order Code: 560-565**

### WEINGART ANGLED

Serrated, tapered tips allow for a firm grip on the wire. Tips are rounded for safety and comfort, yet small enough to fit between brackets. Angled at 45° for access to hard to reach areas. Not for use on hooks or crimpable stops.

**Order Code: 560-008-2**  
 Made in USA  
**Order Code: 500-008**  
 Made in Poland  
**Order Code: 570-008**  
 Made in China

### CRIMPING PLIER

For crimping stops and hooks onto archwires.

**Order Code: 570-563** Straight  
**Order Code: 570-564** Angled

**E. A. Beck & Co • Order Toll Free: (800) 854-0153**

ph (949) 645-4072 • fx (949) 645-4085 • 657 W 19th St #E • Costa Mesa, CA 92627

email: [EABeck@sbcglobal.net](mailto:EABeck@sbcglobal.net) • [www.beckinstruments.com](http://www.beckinstruments.com)

## PLIERS



**UNIVERSAL STEP UP/STEP DOWN PLIER**  
Allows you to make various 1st order and 2nd order (by-pass) bends intra-orally without removing the entire arch wire. One plier makes two .030 offset bends. One tip makes a step up; by reversing the plier, it makes a step down. OTHER SIZES AVAILABLE.

Order Code: 560-575C  
Made in the USA  
Order Code: 560-575p  
Made in Pakistan



**RICKETTS LOOP PLIER**  
A longer length loop and spring bending plier than the standard Tweed loop. Wire sizes for the first turret is .016", the second is .018" and the third is for .022" wire. One beak is finely serrated to prevent wire slippage. This plier may also be used to bend or adjust soft .016" rectangular arch wire.

Order Code: 500-027a  
Made in Poland  
Order Code: 570-027  
Made in Japan



**BRACKET REMOVING PLIERS**  
500-345 straight  
500-344 Angled  
Regular price \$111.00



**POSTERIOR BAND REMOVER**  
Advanced tip design for easier band removal with less force on the tooth. Excellent fulcrum advantage.

Order Code: 500-355



**NiTi THREE PRONG**  
Particularly fine and delicate tips for all thermal and NiTi-arch wires. For hard wire up to .4mm / .016".

Order Code: 510-025



**NiTi DISTAL BENDER**  
Bend NiTi wires intraorally  
Order Code: 500-004

## LIGATURE DIRECTOR & SEATERS

**LIGATURE DIRECTOR "E"**  
For directing and tucking ligatures. Stainless steel.  
Order Code: 200-656



**LIGATURE TUCKER**  
Blacks #1. For tucking ligatures. Stainless steel.  
Order Code: 200-648



**AND SEATER**  
Large pattern, also used for lapping bands. All stainless steel.  
Order Code: 500-513

**BRACKET HEIGHT GAUGE**  
Order Code: .022" Slot 500-310  
.018" Slot 500-311



**BOONE GAUGE**  
For measuring bracket height.  
Order Code: 500-320

# MEASURING DEVICES


**DIGITAL CALIPER**  
 Order Code: 100-975M



- Includes 1.5V button cell battery
- Fluid-resistant stainless steel frame
- Accurate to 0.001"
- 4 way - inside, outside, depth and step measurements
- Precision ground jaws
- Metric and SAE digital display - .001" / 0.1 mm readings
- Can toggle between metric and inch measurement
- On/Off switch with automatic shutoff



**RULER**  
 White 6" flexible ruler with millimeter and inch measurements.  
 Order Code: 800-001




**CEMENT**  
 "Ultra Band-Lok® is intended for use as an orthodontic band cement an for occlusal buildups".  
 Order Code: 800-110

# HEMOSTATS AND SCISSORS




**BECK MATHIEU FOR ELASTICS**  
 Using wire with this plier voids the warranty.  
 Order Code: 500-081

Positive locking pliers with fine serrated beaks for holding and placing elastic materials without tearing them. Imported from Germany. Lifetime warranty against defects in material and workmanship. The Mathieu for Elastics is for elastics only.



**MOSQUITO FOR ELASTICS**  
 These straight latching hemostats are ideal for placing elastic materials and are guaranteed not to tear the elastics nor allow them to slip.  
 Order Code: 500-231



**BECK MATHIEU FOR TWISTING LIGATURES**  
 Order Code: 500-083



**RICKETTS SCISSOR**  
 A dual purpose instrument. Not only an unusually fine ligature and band scissor, but with serrated tips and blade sides. May be used as a tucker as well.  
 Order Code: 500-229

**INVISIBLE RETAINER SCISSORS**  
 Scissors for foils and splints. Curved with serrations. Chrome plated. 95mm / 3.74". May not be sterilized.  
 Order Code: 510-215



**LIGATURE SCISSOR**  
 For trimming invisible retainer material, ligature wire and bands. Please specify serrated, if desired. Made in Germany of surgical stainless steel. Lifetime warranty against defects in material and workmanship.  
 Order Code: Cvd Heavy 500-215

## CUTTERS



### DISTAL END CUTTERS

Cut and hold cutters with carbide inserts.

Order Code: 500-135

Made in Poland



### HARD WIRE CUTTERS

For hard wires up to .020"

Order Code: 500-140

Made in Poland

Order Code: 580-139

Made in USA

### WIRE SIZE LIMITS

- Rectangular .022" x .025"
- Round .020"
- Braided .016"



### LIGATURE CUTTERS

Carbide inserts - For soft ligature wire up to .012".


Order Code: 500-138

Made in Poland

Order Code: 580-138

Made in USA

We guarantee when our cutters are used with the recommended wire limits they will stay sharp for the first year or we will sharpen them free of charge!




### DIAGONAL WIRE CUTTERS

A tested and proven cutter. Guaranteed not to nick. Ideal for laboratory use. For hard wire up to .036". Chrome plated - may not be steam autoclaved.

Order Code: 500-105

## DISTAL BENDER & BRACKET TWEEZERS




### DISTAL BENDER

Order Code: Single ended .022" 650-300

Order Code: Single ended .030" 650-330

Order Code: Double ended 650-360

- Makes cinchbacks a cinch to do
- Available in 2 sizes .022" and .030"
- Works great on all SS wires, TMA and most Nitinol wires
- Bends distal ends of arch wires
- Bends tails of springs



### KEAT MOLAR BRACKET TWEEZER

Tip of Keat fits through an .018" - .022" tube for easy transfer. Amazing control and vision for placing buccal tubes and brackets. Bottom of tapered handle can be used as a slot aligner. Eliminates slippage and shooting. These tweezers have extra delicate tips and are not guaranteed against breakage.

Order Code: Regular Handle 5 3/4"

600-023k



### TWEEZER STYLE BRACKET PLACERS

Order Code:

600-032



### CROSS ACTION BRACKET PLACERS

Order Code: 4 1/2"

319a

Order Code: 5 1/4"

320a

# DRY FIELD & MINI MOLD KITS

## COTTON MOUTH™

Order Code: 800-200

The Expanded Duty Cotton Mouth™ system provides a dry field for numerous dental and orthodontic procedures while affording your patient total comfort. Suggested uses:

- Power Bleaching
- Cosmetic Bonding
- Porcelain Veneers
- Composite Restoration
- Orthodontic Bonding
- Anterior Crown and Bridge
- Lantophonesis



Expanded Duty Cotton Mouth Dry Field System


Kit includes:

- 1 Expanded Duty Bite Block
- 2 Short Aspirator Tips
- 3 Long Aspirator Tips
- 4 Pads 5mm
- 4 Pads 7mm
- 25 Disposable Aspirator Tubes

## MINI MOLD KIT

Order Code: 500-100

Light-cure your own resin attachments in seconds! Make fixed retainers, tongue cribs, and bite-opening appliances without impressions. Mini-Molds do the job quickly and easily with 7 different reusable tips. Just fill the tip (choose from 7 different treatment aids) place on tooth, and cure. In seconds, Mini-Mold tips precisely cast a light cured resin appliance.



Mini-Mold Starter Kit includes:

- |                       |                  |
|-----------------------|------------------|
| 3 Habit Control       | 3 Button Bonder  |
| 3 Anterior Bite-ramp  | 3 Bracket Bonder |
| 3 Posterior Molar Pad | 3 Tube Bonder    |
| 3 Wire Bonder Small   | 2 Bonder Handles |
| 3 Wire Bonder Large   |                  |

Emergency Curing Light  
(sold separately)

Order Code: 800-510

10 second bond with quick cure  
LED light is not harmful to the eyes  
Includes three AAAA batteries for over 10 hours of operation  
Strength of cure is the same as regular curing light



# ELASTIC POWER THREADS AND TUBING



ELASTIC POWER TUBING  
Clear, Hollow, .030", 50 ft. spool

Order Code: 502-332

ELASTIC POWER THREAD  
Grey, Solid, .025", 50 ft. spool

Order Code: 502-331

SLIP-NOT POWER THREAD  
Slip-Not rough surface holds knots more securely. Clear, .030", 25 ft. spool.

Order Code: 502-333



A.J. WILCOCK® AUSTRALIAN

Stainless steel wire. Offers superior hardness and tensile strength for a multitude of applications where exceptional resiliency is desired. NiTi™ available upon request  
Order Code: 800-320 .016" Special - Higher force with minimal breakage.

Order Code: 800-630 .022" Special + - Requires caution when bending to avoid fracture. 25' per spool.


A.J. Wilcock® Preformed Archwires  
With Intermaxillary Circles

Elastics engagement is simple with A.J. Wilcock® looped preformed archwires. Available in **Special Plus grade wire**. Sold in packs of two. Many sizes available.



A.J. Wilcock® Preformed Arches

Preformed Australian Archwires are unsurpassed in resiliency and hardness with incredible levels of stored energy for rapid tooth alignment. Available in round or rectangular, **Special Plus grade wire**. Sold in packs of 25. Many sizes available



# AIR ROTOR STRIPPING

## Super Flexible Diamond Discs



Grit covers 3mm of disk edge



	.3mm	Double Sided	3/4", 19mm Diameter	7/8", 22mm Diameter
	.2mm	Safe Outside	3/4", 19mm Diameter	7/8", 22mm Diameter
	.2mm	Safe Inside	3/4", 19mm Diameter	7/8", 22mm Diameter

Permanently mounted on mandrels  
For separating and contouring interproximal areas

Double sided: Diamond coating on both sides  
Safe outside: Diamond coating on inside surface  
Safe inside: Diamond coating on outside surface



Diamond Discs with holes

	.15mm	Double Sided	3/4", 19mm Diameter	7/8", 22mm Diameter
	.1mm	Safe Outside	3/4", 19mm Diameter	7/8", 22mm Diameter
	.1mm	Safe Inside	3/4", 19mm Diameter	7/8", 22mm Diameter

## SEPARATING STRIPS

Genuine carborundum galvanically coated on stainless steel. Can be cut to desired length before use. Each strip is 5 7/8". Sold in packages of twelve.



### SINGLE SIDED

- 102-001 Fine 2mm
- 102-002 Fine 4mm
- 102-005 Medium 2mm
- 102-006 Medium 4mm
- 102-007 Medium 6mm
- 102-008 Medium 8mm

### SINGLE SIDED

- 102-011 Coarse 2mm
- 102-012 Coarse 4mm
- 102-013 Coarse 6mm
- 102-014 Coarse 8mm

### DOUBLE SIDED

- 102-015 Medium 6mm
- 102-016 Medium 8mm

# RETRACTORS & MARKERS



**DOUBLE ENDED LIP RETRACTORS**  
Features both a pedo and teen/adult size. Comfortable and easy to hold, yet allows an excellent field of vision for treatment and photographs. Can withstand temperatures of 370 F. Amber tint identifies it as an autoclavable plastic instrument.

Order Code: 600-027



## ARCH MARKER DISPENSER

Order Code: 800-445



## DISPOSABLE ARCH MARKERS

Order Code: Red 800-446  
White 800-447



## HIGH HEAT BITE STICKS

For high heat up to 440 F. Ser-rated stainless steel insert ensures positive contact without slipping. Sold individually.

Order Code: 800-100



## MAGENTIC BUR HOLDER

Nickel-plated steel. Diameter is 1 7/16" and is 1" in height.

Order Code: 500-013